

# DENV E gene Amplification for Genotyping

(Version.1.2 20240111)

## MATERIALS AND CONSUMABLES

• TransStart®FastPfu Fly DNA Polymerase Kit	TransGen	AP231
• RNase-free H <sub>2</sub> O	-	-
• DL 2,000 DNA marker	Takara	3427
• Primer:		
D1F1: CACATGCCATAGGAACATCCA	D1R1: ATGAGCCTGTGCACATCACA	
D1F2: GGAACAGACAAGATTTGCTGGT	D1R2: TGCCACTTCCACATTTGAGT	
D2F1: TGGCATAACCCATAGGAACGA	D2R1: CCTTTGAGCTGTAGTTTGTCCA	
D2F2: CCTCGACTTCAATGAGATGGT	D2R2: TTGAAGGGGATTCTGGTTGGA	
D3F1: AGGGTTCACAATACTAGCCCTA	D3R1: GGGCTACAACAGAAACACCA	
D3F2: GTTCTCCATTCTGGTTGTCTGA	D3R2: GTTCTCCATTCTGGTTGTCTGA	
D4F1: AGCTGGATACTCAGAAACCCAGGATT	D4R1: ACATCCTGTCTCTTGGCATGAGG	
D4F2: CCTCATGCCAAGAGACAGGATGT	D4R2: AATTTGTACTGTTCTGTCCAAGTGTG	

## PROCEDURES

### I. Total Viral RNA extraction & Reverse transcription

### II. PCR

☐ 1. Prepare PCR reaction as follow,

COMPONENT	AMOUNT (μl)	SAMPLE NUM	OPERATE
Template	2.00	-	-
<input type="checkbox"/> Nuclease-free water	31.00	× _____	
<input type="checkbox"/> Forward primers (10 μmol/l)	1.00		
<input type="checkbox"/> Reverse primers (10 μmol/l)	1.00		
<input type="checkbox"/> 5 × TransStart FastPfu Fly Buffer	10.00		
<input type="checkbox"/> tNTPs (2.5 mmol/l)	4.00		
<input type="checkbox"/> TransStart FastPfu Fly DNA Polymerase	1.00		
TOTAL		50.0	

☐ 2. Place in a thermocycler and run the following program,

STEP	TEMPERATURE	TIME	CYCLE
Pre-denature	98	2 min	1
Denature	98	10 s	40
Anneal	59	20 s	
Extend	72	15 s	
Re-Extend	72	5 min	1
Stop	16	-	-

☐ 3. Analysis the product by Agarose gel electrophoresis.

### III. The PCR products Sanger sequencing